



Bouliga Description

Part One Product Brief

Bouliga brand is made of cross-linked sodium hyaluronate of non-animal origin, a substance naturally present in our body, it is safe and completely biocompatible, as well as wholly biodegradable.

Bouliga has been designed from biofermentation hyaluronic acid, a substance naturally present in our body, it is safe and completely biocompatible, as well as wholly biodegradable.

Function is to treat all sorts of furrow, forehead horizontal grain, eyebrow between vertical grain, dimple grain, nasolabial furrow grain, labial furrow, nasolabial root furrow, cheek furrow, ear before furrow, neck crosswise grain and so on can handle easily.

Part Two Bouliga Features

Sodium hyaluronate is a natural polysaccharide in human skin, which has the characteristics of maintaining skin moisture, nourishing skin and anti-aging. The effective ingredients in sodium hyaluronate fillers are the same as the natural unmodified hyaluronic acid in human body, which can be completely degraded in human body without residue and long-term side effects. Sodium hyaluronate filling agent has been called "God clay" for a long time. It has double functions of filling and shaping. It can fill depressions, smooth wrinkles, moisten and moisturize skin, and the effect is quick and magical. It is a non-toxic, non repellent, non allergic and inflammatory skin filling product.

Bouliga is a non animal source sodium hyaluronate filling agent. It has excellent cosmetic and plastic effect, simple and safe treatment process, excellent compatibility with human body, stable, natural, not removed after injection, and can be absorbed quickly and effectively. Injection is used to achieve the effect of cosmetic surgery and avoid the risk of surgery. Healthy and natural beauty filler products can not only fill static wrinkles and dynamic wrinkles, but also make the skin moisturize for a long time.

Part Three Bouliga Advantage

① Production of hyaluronic acid fillers from high purity hyaluronic acid

The raw material is hyaluronic acid from non animal sources produced by microbial fermentation, and the quality of the raw material is higher than the requirements of European Pharmacopoeia and Japanese Pharmacopoeia. The content of nucleic acid, protein and endotoxin is far lower than the standards of European Pharmacopoeia and Japanese Pharmacopoeia.

② Perfect quality management system

Based on years of experience in quality management system, the company has established an international medical device quality management system. Scientific and perfect process design and process validation to ensure product safety and effectiveness. Strict process control to improve product quality.

③ Lasting effect

Bouliga's unique polymer hyaluronic acid strengthening and stabilizing technology enables better elasticity and viscosity, ensures smoothness and fixity, and is more conducive to support and shaping.

④ High security

Bouliga's production adopts the international strict way to control the quality. Its spatial network structure is more compact and its injection performance is high. Its osmotic pressure is basically the same as tissue osmotic pressure, reducing the incidence of edema.

⑤ Stable

Bouliga's g 'value is close to the rheological characteristics of normal human dermis and subcutaneous tissue. After injection, it has a stable shape and no displacement.

Part Four Bouliga Effect

1. Face and body shaping:

Hyaluronic acid can be powerful in shaping, with a strong nose, sexy lips and perfect face shape To upgrade one's appearance, you can achieve a "gold standard" face shape with just one shot, without surgery.

Beauty type range: nose augmentation, lip augmentation, chest augmentation, buttock augmentation, temple augmentation, cheek augmentation, earlobe augmentation, forehead augmentation, whole face adjustment, eyebrow augmentation, chin augmentation

2. Strong wrinkle removal:

Hyaluronic acid can smooth and eliminate large wrinkles and small wrinkles. After the injection of hyaluronic acid, it can fill the concave holes and wrinkles in the cortex, restore the skin surface to a full state, instantly smooth the wrinkles, and make people instantly have the natural beauty of youth.

Wrinkle removal scope: forehead transverse lines, eyebrow vertical lines, dimple lines, nasolabial groove lines, lip wrinkles, nasal root transverse lines, cheek wrinkles, ear front wrinkles, neck transverse lines It's all easy.

3. Water locking and moisturizing:

Hyaluronic acid is abundant in deep skin, and the content of hyaluronic acid is the highest in infant skin. 1g of hyaluronic acid can absorb 500cc of water, with 500 times of water absorption capacity, making the skin no longer "thirsty". After absorbing water, hyaluronic acid will make the elastic fiber and collagen in the skin stay in a moist environment for a long time, and instantly restore the skin's elasticity and moisturize the skin.

Part Five Specifications

1ml 2ml 5ml 10ml 20ml 50ml

Part Six Injection Method

Linear injection, multi-point sequence injection, fan-shaped injection, staggered injection

Part Seven Attention

- ① Before injection, doctors plan treatment plan, remove makeup and clean the patient's face, apply ice, disinfect the face, apply anesthetic ointment and inject.
- ② Prohibited for patients with a history of severe allergies or autoimmune diseases
- ③ It is forbidden to inject into existing permanent fillers or different kinds of fillers
- ④ It is forbidden to use in the parts with inflammation or skin diseases, such as infection, psoriasis and herpes.
- ⑤ It is forbidden to use in patients with abnormal coagulation function, or patients who have received anticoagulant treatment (such as aspirin) within 5 days, and who are taking blood activating and stasis removing drugs.
- ⑥ After opening the package, please use it all at one time. Do not use it again next time if you haven't used it up . Do not share the needle with others.
- ⑦ Do not touch the injection area within 6 hours after injection, do not use cosmetics and other irritant articles within 2 hours after injection, do not press hard 48 hours after injection, avoid large facial movements such as laughing and crying, and do not touch the injection site with water or local pollution.
- ⑧ Avoid exposure to high temperature, sauna or extreme cold before swelling and redness subside.
- ⑨ Avoid spicy and irritant food and alcohol within one week after injection.
- ⑩ Injection is prohibited during menstruation, lactation and pregnancy.

Part Eight Shape

The product is a viscous suspension containing transparent or translucent particles, which is composed of cross-linked sodium hyaluronate from non animal sources and a small amount of free sodium hyaluronate. It is packed in a disposable glass syringe sealed by a cone head protective cap. The syringe interface is a rotary lock interface.